

Product Data Sheet

Product Code: CLN4T

Issued: September 28, 2016

What is CleanSave 4T?

A powerful detergent for toilets and urinals. It cleans and brightens porcelain, ceramic and stainless steel sinks, basins, urinals and toilet bowls.

CLEANSAVE 4T is a part of the CLEANSE Closed Loop Chemical and Dispensing System which prevents cleaners from having contact with the concentrated chemicals.

CleanSave Range

CLEANSAVE is a range of performance driven, concentrated chemicals that are supplied in tamper-proof packaging for dispensing automatic dilution systems. Each convenient 2.4L bottle comes fitted with a safety insert that stops cleaners from being able to access the concentrate. The products can only be diluted through dispensers which provides complete control over dilution, stopping over-dosing and product wastage and improving safety.

Key Benefits

- ✓ Highly concentrated to get more out of each bottle
- ✓ Provides control over dilution
- ✓ Virtually no product wastage or 'over dosing'
- ✓ Excellent day-to-day toilet cleaning performance
- ✓ Rapidly removes soilage
- ✓ Colour and number coded for ease of use and safety
- ✓ Environmentally responsible, biodegradable and phosphate-free
- ✓ Reduces packaging, transport and storage requirements

Environmental Care

CLEANSAVE 4T has been made from ingredients that have been chosen to give the following properties:



- ✓ Highly biodegradable surfactants
- ✓ pH mildly acidic
- ✓ Phosphate-free
- ✓ Non-flammable
- ✓ Non-toxic
- ✓ Naturally derived surfactants
- ✓ Minimal VOC's (volatile organic compounds)
- ✓ Recyclable packaging

How Does It Work?

CLEANSAVE 4T's special formula will dissolve rust but will not affect stainless steel, even during periods of prolonged contact.

For Use On...

CLEANSAVE 4T is a mild-acid cleaner for toilets, urinals, tiled washroom floors, walls, baths, basins and shower recesses. It dissolves all types of deposits from porcelain, ceramic and stainless steel surfaces, and leaves surfaces perfumed with a eucalyptus fragrance.

Caution

- CLEANSE 4T is an acid-based cleaner. Avoid contact with skin and eyes.
- Because it is an acidic product, CLEANSE 4T is not suitable for use on marble, terrazzo, Formica or Laminex surfaces. (Even a slice of lemon can damage these surfaces.)
- NEVER MIX CLEANSE 4T with other chemicals. Never mix with BLEACH.

Technical Data

Properties



COLOUR – Transparent blue liquid
ODOUR – Eucalyptus fragrance
pH = Undiluted pH is 2-3
FOAM – Thick foam with good foam stability

Health

CLEANSAVE 4T is not classified as Dangerous Goods according to the criteria of the Australian Safety and Compensation Council (ASCC). It contains no ingredients that are classed as Toxic.

Quality

The design, manufacture and supply of all Agar chemical products is controlled by the Agar Quality Management System which is registered and externally audited by SAI Global as complying with the requirements of AS/NZS ISO 9001 "Quality Management Systems – Requirements".

First Certified: 30 April 1996

SAI Certificate No.: QEC7358

Application

The ready-to-use solution of CLEANSE 4T is available from your Agar Detergent Dispenser for bottle fill. Simply select the product and push the button while holding the hose in your bottle.

Pour 60mls of the solution around the rim of the toilet bowl and allow it to run down, covering the whole bowl. Scrub with a nylon brush then flush after a few minutes.

Agar Cleaning Systems maintains Safety Data Sheets (SDS) on all of its products. These sheets contain information that you may need to protect your employees and customers against health or safety hazards associated with our product. Agar Cleaning Systems recommends that you obtain a copy of the respective SDS sheet prior to using this product. The information in the Product Data Sheet is based on data we believe to be reliable. It is offered in good faith, but without guarantee, as conditions and methods of use of our product are beyond our control.